


T2J

**TONGUE, LIPS & JAW
APPLIANCE**

**STRENGTHENS
TONGUE MUSCLES**

**STRENGTHENS
THE LIP MUSCLES**

**IMPROVES
LIP SEAL**



**ASSISTS IN
ESTABLISHING
NASAL
BREATHING**

myoTALEA[®]
TONGUE AND LIP
EXERCISE APPLIANCE


UNDERSTANDING THE PROBLEM

Patients with orthodontic, jaw, breathing and sleep problems lack sufficient strength in the tongue, lip and jaw muscles. These muscles need to be strengthened to resolve the underlying causes of sleep-related breathing problems.

The *Myotalea*® *TLJ* appliance is an active myofunctional appliance used for strengthening the tongue, lip and jaw muscles. Additional exercises also target the throat and airway muscles, which are typically weak in patients who show symptoms of Breathing Disordered Sleep (BDS). This is an essential part of myofunctional orthodontic, airway and breathing treatments.


The *TLJ* has numerous features that will aid in the correction of dysfunctional soft tissue habits including the following:

- Habitual mouth breathing
- Lips apart rest posture
- Low tongue rest posture
- Incorrect swallowing pattern
- Lack of strength in throat and airway muscles



INDICATIONS FOR USE

The *Myotalea*® TLJ appliance provides treatment for patients with symptoms of sleep-related breathing disorders and jaw joint dysfunction. Patients who present poor muscle tone in the tongue, lip, jaw, throat or airway areas are suitable for *Myotalea*® treatment.



FLOW SEQUENCE

If advised by your practitioner, these exercises can be combined to form an exercise flow sequence as follows:

Simultaneous
Tongue Press & Hold
+ Lip Press & Hold
(3 breaths)



Tongue Press & Hold
(3 sets)




Lip Press & Hold
(3 sets)



Jaw Press & Hold
(3 sets)




Simultaneous
Tongue Press & Hold
+ Lip Press & Hold
(3 breaths)



ADVANCED TECHNIQUE: **HEAD TILT VARIATION**

Repeat the steps with the head tilt posture. There should be no discomfort in head, neck and jaws during this variation.



TREATMENT WITH MYOTALEA® TLJ


Malocclusion, sleep-related breathing problems and jaw joint disorders are all indications of poor orofacial muscle strength. These tongue, lip, jaw, throat and airway muscles are typically underused during growth years, which subsequently affects growth and development.

Patients who suffer from Sleep Disordered Breathing or Temporomandibular Joint (TMJ) Disorder are often mouth breathers. Habitual mouth breathing and poor muscle tone leads to airway collapse while the patient sleeps and must be addressed to achieve long-term symptom relief.

The *Myotalea*® TLJ is to be used as an aid in alleviating Sleep Disordered Breathing problems.

DESIGN FEATURES

- 1 **Lip Press Tube** strengthens the lip muscles, which improves lip seal and nasal breathing.
- 2 **Tongue Press Tube** improves tongue strength and position.
- 3 **Air Springs** improve function of the jaw joint and muscles.



TONGUE, LIPS & JAW
APPLIANCE



INSTRUCTIONS FOR USE

It is recommended that the *TLJ* appliance is initially used for approximately **3 minutes twice a day** to gain the full benefits of its exercising features. This simple sequence will aid in developing better muscle strength and tone. Build up to the recommended time durations based on your current functional abilities. It is important that you challenge yourself and push your muscles to fatigue to start building strength and tone.

TONGUE

Tongue Press & Hold

- Place the *TLJ* appliance in the mouth with the Tongue Press Tube facing upwards.
- Compress the Tongue Press Tube between the tongue and the roof of the mouth.
- Hold for 3 breaths while keeping the lips together and breathing lightly through the nose.
- Head tilt variation*

LIPS







Lip Press & Hold

- Place the *TLJ* appliance in the mouth with the Lip Press Tube between the lips.
- Press lips firmly together and compress the Lip Press Tube in the process.
- Hold for 3 breaths while keeping the lips together and breathing lightly through the nose.
- Head tilt variation*

JAW

Jaw Press & Hold

- Place the *TLJ* appliance in the mouth with the teeth comfortably placed.
- Lightly bite down on both sides of the appliance to compress the Air Springs.
- Hold for 3 breaths while keeping the lips together and breathing lightly through the nose.
- Head tilt variation*



***HEAD TILT VARIATION:** Tilt the head backwards while performing exercises to improve strength in the throat and airway muscles.

DESIGN FEATURES

- **Lip Press Tube**
strengthens the lip muscles, which improves lip seal and encourages nasal breathing.
- **Tongue Press Tube**
improves tongue strength and position through elevation and compression techniques.
- **Air Springs**
improve the function of the jaw joints and muscles.
- **Anterior Lip Bumpers**
discourage overactive muscles in the chin and lip area and can assist in establishing correct swallowing habits.
- **Tongue Elevators**
encourage the tongue to remain in the correct elevated position.
- **The Myotalea® TLJ appliance**
is effective for both isometric and isotonic exercising of the orofacial musculature.


CLEANING THE TLJ: Clean your *Myotalea® TLJ* appliance under warm running water every time you remove it from your mouth. Use *Myoclean™* tablets to correctly clean once every week. Ask your Doctor about *Myoclean™*, the recommended cleaning agent for all *MRC* appliances.

CAUTION: CHILDREN MUST BE SUPERVISED BY AN ADULT WHEN USING THE TLJ.

Lifetime use: 6 months based on everyday wear.
Replace after 12 months regardless of wear time.

 **Myotalea® by Myofunctional Research Co.**
Myofunctional Research Co. Australia:
44 Siganto Drive, Helensvale QLD 4212 Australia
Tel: +61 7 5573 5999

USA Representative:
Myofunctional Research Co. USA
9267 Charles Smith Avenue, Rancho Cucamonga CA 91730 USA
PO Box 2781, Rancho Cucamonga CA 91730 USA
Tel: +1 866 550 4696

 **European Representative:**
Myofunctional Research Europe B.V.
Gompensstraat 21c 5145 RM Waalwijk The Netherlands
P.O. Box 718, 5140 AS Waalwijk The Netherlands
Tel: +31 416 651 696

Single Patient Use.



**MYOFUNCTIONAL
RESEARCH CO.**

**INNOVATIVE DENTAL APPLIANCE
TECHNOLOGY AND EDUCATION**

